SISTEMAS DE IMPLANTES IMPLANT SYSTEMS

www.gmidental.com

frontier

frontier
PEAK

Implantes dentales de conexión interna hexagonal
Hexagonal internal connection dental implants

Made in Spain

Índice - Summary

CARACTERÍSTICAS GENERALES

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- Kit quirúrgico **gmi frontier/avantgard**
- Kit quirúrgico **gmi PEAK**

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KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

IMPLANT SYSTEMS

- **gmi frontier** dental implants 8
- Insertion sequence for **gmi frontier** implants 9
- **gmi frontier** implants range 10
- Drilling sequences for **gmi frontier** implants 12
- **gmi frontier PEAK** dental implants 14
- Insertion sequence for **gmi frontier PEAK** implants 15
- **gmi frontier PEAK** implants range 16
- Drilling sequences for **gmi frontier PEAK** implants 17

ATTACHMENTS

- **gmi frontier** healing abutments 18
- **gmi frontier** millable-implant abutments 19
- **gmi frontier** esthetic angled abutments 20
- **gmi frontier** straight abutments 20
- **gmi frontier** temporary cylinders 21
- **gmi frontier** titanium base 22
- **gmi frontier** ball abutments 23
- **gmi frontier** multi-esthetic abutments 24
- **gmi frontier** bioesthetic bases 27
- **gmi frontier** replicas 29
- **gmi frontier** impression coping 30
- **gmi frontier** castables 31
- **gmi frontier** CrCo bases 32
- **gmi frontier** clinic and laboratory screws 32
- **gmi frontier** MEA and BE screws 33
- Recommended tightening torques 34

GMI SURGICAL KITS

- **gmi frontier/avantgard** surgical kit 37
- **gmi PEAK** surgical kit 38

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

To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.





SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

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The **gmi** **frontier** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

		
ø3,30-ø3,75 mm	ø4,25 mm	ø4,75-ø5,75 mm
Incisivos (Excepto centrales maxilares) Incisors (Except central maxillary incisors)	Caninos • Canines Premolares • Premolars Incisivos centrales maxilares Central maxillary incisors	Molares Molars

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The **gmi** **frontier** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi** **frontier** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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

To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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In the chart are detailed all the symbols appearing on the **gmi frontier** implant system's labelling and packaging and their corresponding description.


Caja - Packaging
Connection and platform

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* Minimum requirements for downloading and opening the electronic instructions for use:

- Stable internet connectivity.
- A device with a modern operating system that supports PDF viewing.
- Adequate storage space on your device.
- A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.
- A PDF reader application installed on your device.

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
gmi frontier internal connection hexagonal implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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① Platforms
RP: Standard $\varnothing 3.30$ mm
WP: Wide $\varnothing 4.30$ mm

② Length

③ Self threading millings
Variable length
3 millings at 120°
Atraumatic tip

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gmi frontier internal connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi frontier** implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single hex-2.45 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.*

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* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the HEX-2.45 mm implant carrier wrench until a slight retention is noticed, addressing the hex of the key with the faces of the implant hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the crestal or subcrestal level. Remove the key.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove key and suture the incision.

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

The following is a description of the drilling sequences for the different internal connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped drill: 1200 - 1500 rpm.
- Pilot drill: 700 - 900 rpm.
- Step drills:
 - Ø3.00 & Ø3.50mm → 400 - 700 rpm.
 - Ø4.00 & Ø5.10mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.

For drilling, external cooling with saline solution is required.


frontier Ø3,30		
Descripción - Description		Ref.
1 Fresa lanza	Lance-shaped drill	KYFoC1225
2 Fresa piloto Ø2.50	Pilot drill Ø2.50	KYFoC2250
3 Fresa cortical Ø3,65	Cortical drill Ø3.65	KYFoC4336
4 Fresa escalonada Ø2,50/Ø3,00	Step drill Ø2.50/Ø3.00	KYFoC5130


frontier Ø3,75		
Descripción - Description		Ref.
1 Fresa lanza	Lance-shaped drill	KYFoC1225
2 Fresa piloto Ø2,50	Pilot drill Ø2.50	KYFoC2250
3 Fresa escalonada Ø2,50/Ø3,00	Step drill Ø2.50/Ø3.00	KYFoC5130
4 Fresa cortical Ø4,00	Cortical drill Ø4.00	KYFoC4340
5 Fresa escalonada Ø3,00/Ø3,50	Step drill Ø3.00/Ø3.50	KYFoC5135

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

gmi frontier PEAK internal connection hexagonal implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the conical core and progressive thread profile facilitate implant insertion and act as a bone compactor allowing high primary stability to be obtained in sparse bones.

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② Cortical thread
1.50 mm

③ Length

④ Self threading millings
Variable length
2 millings at 180°
Atraumatic tip
Cut in both directions

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gmi frontier PEAK internal connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi frontier PEAK** implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single hex-2.45 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.*

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

The following is a description of the drilling sequences for the different internal connection implant models, as well as the recommended drilling conditions for their use:

- Ball and initial drill: 1000 rpm.
- Step and conical drills:
 - Ø2.00 / Ø2.80 mm → 500 - 700 rpm.
 - Ø3.20 / Ø3.65 mm → 500 - 700 rpm.
 - Ø3.70 / Ø4.15 mm → 400 - 700 rpm.
 - Ø4.15 / Ø4.65 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.


For drilling, external cooling with saline solution is required.

frontier PEAK Ø3,75		
Descripción - Description		Ref.
1 Fresa inicial	Initial drill	KYFoC1401
2 Fresa escalonada Ø2,00/Ø2,80 mm	Step drill Ø2.00/Ø2.80 mm	KYFoC5128
3 Fresa cónica Ø3,20 mm	Conical drill Ø3.20 mm	KYFoC7032
4 Fresa cortical Ø3,65 mm	Cortical drill Ø3.65	KYFoC4337
5 Fresa escalonada Ø3,20/Ø3,65 mm	Step drill Ø3.20/Ø3.65 mm	KYFoC5137

frontier PEAK Ø4,25		
Descripción - Description		Ref.
1 Fresa inicial	Initial drill	KYFoC1401
2 Fresa escalonada Ø2,00/Ø2,80 mm	Step drill Ø2.00/Ø2.80 mm	KYFoC5128
3 Fresa cónica Ø3,20 mm	Conical drill Ø3.20 mm	KYFoC7032
4 Fresa cónica Ø3,70 mm	Conical drill Ø3.70 mm	KYFoC7037
5 Fresa cortical Ø4,15 mm	Cortical drill Ø4.15 mm	KYFoC4342
6 Fresa escalonada Ø3,70/Ø4,15 mm	Step drill Ø3.70/Ø4.15 mm	KYFoC5142

frontier PEAK Ø4,75		
Descripción - Description		Ref.
1 Fresa inicial	Initial drill	KYFoC1401
2 Fresa escalonada Ø2,00/Ø2,80 mm	Step drill Ø2.00/Ø2.80 mm	KYFoC5128
3 Fresa cónica Ø3,20 mm	Conical drill Ø3.20 mm	KYFoC7032
4 Fresa cónica Ø3,70 mm	Conical drill Ø3.70 mm	KYFoC7037
5 Fresa cónica Ø4,20 mm	Conical drill Ø4.20 mm	KYFoC7042
6 Fresa cortical Ø4,65 mm	Cortical drill Ø4.65 mm	KYFoC4347
7 Fresa escalonada Ø4,15/Ø4,65 mm	Step drill Ø4.15/Ø4.65 mm	KYFoC5147

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
Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Alternative tightening with a flat tip screwdriver.
- Available in emergency heights from 2.00 to 6.00 mm.
- Apical bit to facilitate the start of threading.

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
When the use of a custom abutment is required, millable abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, millable implant abutments can be used like impression coping, threading them to the implant screw through the short impression screw and using them to make the negative of the patient's mouth for subsequent positivation on plaster.

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Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 4.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Clinic screw ref.: KDTOC3601 included.



Plataforma Platform	Ref.
RP ●●	KDBoF3603
WP ●	KDBoF3602

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw ref: KDT0C3601 included.
- Abutment shoulder adapted to the shape of the gum.

Pilares inclinados estéticos Esthetic angled abutments		
Plataforma Platform	Ang. (°)	Ref.
RP ●●	15°	KDMoF3634
	20°	KDMoF3635
WP ●	15°	KDMoF3636
	20°	KDMoF3637

Pilares inclinados Angled abutments		
Plataforma Platform	Ang. (°)	Ref.
RP ●●	20°	KDMoF3613
WP ●	20°	KDMoF3614

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Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Abutments in b/f 2.45 mm hex non-rotary version.
- Clinic screw ref: KDT0C3601 included.
- Available in three emergency heights from 1.50 to 3.50 mm.

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When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary cylinder is used, defined as technical plastic machined attachment that once attached to the crown of the implant with clinic screw, will absorb most of the masticatory stresses while avoiding overloading the implant while the osseointegration process is completed.

Key features

- Machined in high strength technical plastic.
- Abutment shoulder adapted to the shape of the gum.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw ref.: KDTOC3601 included.


Plataforma Platform	Tipo conexión Connection type	Ref.
RP ● ●	Anti-rotatoria • Non-rotary	KDHof3605
	Rotatoria • Rotary	KDHof3607
WP ●	Anti-rotatoria • Non-rotary	KDHof3606
	Rotatoria • Rotary	KDHof3608

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Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw ref.: KDT0C3601 included.
- Available in emergency heights of 0.60 mm, 1.50 mm and 3.00 mm.
- Non-rotary titanium base with indexation.
- Scan directly on Titanium base or alternative scan on Scan abutment.

Base titanio - Titanium base			
Plataforma Platform	Altura (mm) Height (mm)	Tipo conexión Connection type	Ref.
RP ●●	0,65*	Anti-rotatoria • Non-rotary	KDMoF3615 ①
	0,65*	Rotatoria • Rotary	KDMoF3616 ①
	1,50	Anti-rotatoria • Non-rotary	KDMoF3638 ②
	1,50	Rotatoria • Rotary	KDMoF3640 ②
	3,00	Anti-rotatoria • Non-rotary	KDMoF3639 ②
	3,00	Rotatoria • Rotary	KDMoF3641 ②
WP ●	0,65*	Anti-rotatoria • Non-rotary	KDMoF3617 ①
	0,65*	Rotatoria • Rotary	KDMoF3618 ①
	1,50	Anti-rotatoria • Non-rotary	KDMoF3642 ②
	1,50	Rotatoria • Rotary	KDMoF3644 ②
	3,00	Anti-rotatoria • Non-rotary	KDMoF3643 ②
	3,00	Rotatoria • Rotary	KDMoF3645 ②

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
* Use the Healing abutment h=3 mm RP (KDPOC3609) and WP (KDPOC3613)

Calcinable - Castable			
Plataforma Platform	Altura total (mm) Total Height (mm)	Cantidad (Ud.) Quantity (Un.)	Ref.
RP ●●	11,00	10	KDCoF3605
WP ●	11,00	10	KDCoF3606

Base titanio gmi frontier - gmi frontier titanium base

Pilar escaneado - Scanabutment

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	5,30	KDIOF3601
WP ●	5,30	KDIOF3602

Cuerpo pilar escaneado - Scanbody

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP - WP ●●	7,50	KDIOC4003

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
Once the healing abutment has formed the mucosal tunnel and the healing abutment has been removed from this location, the corresponding ball abutment is used, which is defined as a machined attachment being fixed to the top of the implant and serving as the retention element for the removable prosthesis, with the previous addition of the retention caps that will be fixed to the top spherical part of the ball abutment.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- Pink retention cap and metal container included.
- Protection disk included.
- Tightening with b/f 2.45 mm hex wrench, part no.: KYLOC0065.
- Instruments for caps: Insertion / removal tool for caps (RH485IC).

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	1,00	KDFoF3609
	2,00	KDFoF3610
	3,00	KDFoF3611
	4,00	KDFoF3612

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.50 mm.
- One-piece straight abutments in rotary model.
- Tightened with hex keys:
RP: B/F 2.00 mm (KYLoCo149)
WP: B/F 2.50 mm (KYLoCo098)
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium as per ISO-5832 with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.

Pilares multi-estético gmi frontier - gmi frontier multi-esthetic abutments


Pilares multi-estético inclinados Multi-esthetic angled abutments

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP ●●	17°	2,50	KDGoF3609
	17°	3,50	KDGoF3610
	30°	4,00	KDGoF3611

Pilares multi-estético rectos Multi-esthetic straight abutments

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	1,00	KDGoC3601 ①
	1,50	KDGoC3620 ②
	2,50	KDGoC3621 ②
	3,50	KDGoC3622 ②
	4,50	KDGoC3623 ②
WP ●	1,00	KDGoC3605 ①
	1,50	KDGoC3624 ②
	2,50	KDGoC3625 ②
	3,50	KDGoC3626 ②
	4,50	KDGoC3627 ②

Pilares cicatrización PME Healing abutments MEA

Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP ●●	Ø 4,80	KDPoF0035
	Ø 6,00	KDPoF0037
WP ●	Ø 6,00	KDPoF0036

Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

Plataforma Platform	Ref.
RP ●●	KDSoF0015
WP ●	KDSoF0016

Pilares multi-estético gmi frontier - gmi frontier multi-esthetic abutments

(x5)

Toma impresión cubeta cerrada (CC) PME Closed tray (CT) impression coping MEA

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
Once the insertion of the implant in the bone bed has been completed, the bioesthetic base is used, which is defined as the machined attachment that, once fixed to the coronal part of the implant, allows single restorations at the abutment level.

Placing the bioesthetic base during surgery allows the soft tissues to heal around it.

Once the tissue has healed, it serves as the basis for the rehabilitation, avoiding the connection and disconnection of the components to the implant and, therefore, keeping the biological seal intact.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- Straight two-piece bases, including the screw.
- Tightened with hex wrench: B/F 2.00 mm (KYLoCo149).
- Impression coping with screw included.
- Castable.
- Tightening torques: BE base (30 N-cm) and healing abutment (Manual: Max. 15 N-cm).

Bases bioestéticas Bioesthetic bases		
Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	1,00	KDGoF3612
	2,00	KDGoF3613
	3,00	KDGoF3614
	4,00	KDGoF3615
WP ●	1,00	KDGoF3616
	2,00	KDGoF3617
	3,00	KDGoF3618
	4,00	KDGoF3619

Pilar cicatrización BE BE healing abutment		
Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP - WP ●●	Ø 4,80	KDPoC0038

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




Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.

Réplicas de implante Implant replicas

Plataforma Platform	Ref.
RP ●●	KDRoC3601
WP ●	KDRoC3602

Réplicas digitales Digital replicas

Plataforma Platform	Ref.
RP ●●	KDRoF3603
WP ●	KDRoF3604

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Plataforma Platform	Ref.
RP ●●	KDRoC0012
WP ●	KDRoC0013

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Plataforma Platform	Ref.
RP ●●	KDRoF0014
WP ●	KDRoF0015

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
Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, are used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- CT impression coping whole formed by CT impression coping + hexagonal abutment + screw.

Tomas impresión cubeta abierta (CA) Open tray (OT) impression coping		
Plataforma Platform	Tornillo Screw	Ref.
RP ●●	Corto • Short	KDSOf3608
	Largo • Long	KDSOf3609
WP ●	Corto • Short	KDSOf3610
	Largo • Long	KDSOf3611

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Tomas impresión cubeta cerrada (CC) Closed tray (CT) impression coping

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
To replace a dental crown with a restoration dental prosthesis, a castable is used, which is defined as a meltable plastic machined attachment that, once positioned over the coronal part of the internal connection implant, will be the initial basis to model the form of the final structure, to be cast in the corresponding metal.

Key features

- Machined castables in white meltable plastic material.
- Physiologically inert and excellent dimensional stability.
- Available in rotary and non-rotary version.
- Clinic screw ref.: KDT0C3601 included.

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
The CrCo bases consist of a machined CrCo interface and a castable chimney. This configuration allows for greater dimensional stability at the contact surfaces of the implant-abutment connection.

They are indicated for single or multiple screw-retained restorations either by conventional cast-on technique or by laser welding.

Plataforma Platform	Tipo de conexión Connection type	Ref.
RP ●●	Rotatoria • Rotary	TS30017
	Anti-rotatoria • Non-rotary	TS31017
WP ●	Rotatoria • Rotary	TS30018
	Anti-rotatoria • Non-rotary	TS31018

Tornillos clínica y laboratorio gmi frontier - gmi frontier clinic and laboratory screws

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
Clinic screws are defined as machined devices used for the final anchoring of the abutments, the millable implant abutment or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of castables to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo clínica Clinic screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ●●	○ 1,20 mm	KDToC3601
	⊗ TX	KDToC3609

Tornillo laboratorio Laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ●●	○ 1,20 mm	KDToC3602
	⊗ TX	KDToC3610

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MEA and BE clinical screws are defined as machined devices used for the final anchoring of the ME and BE abutments, while laboratory screws are the machined devices used for the temporary anchoring of ME and BE abutments to the replica during the fabrication of the prostheses in the laboratory.

Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm and 2.00 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo PME MEA screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP ●●	⬡ 1,20 mm	KDToC0056
	⬡ TX	KDToC0070
WP ●	⬡ 1,20 mm	KDToC0060
	⬡ TX	KDToC0071

Tornillo laboratorio PME MEA laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP ●●	⬡ 1,20 mm	KDToC0057
	⬡ TX	KDToC0072
WP ●	⬡ 1,20 mm	KDToC0058
	⬡ TX	KDToC0073

Tornillo BE BE screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ●●	⬡ 1,20 mm	KDToC4016
	⬡ TX	KDToC0074

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
gmi frontier/avantgard surgical kits consist of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi frontier** and **avantgard** dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134 °C.

To complement this kit, **gmi** has the **PEAK drill Kit**, allowing it to be easily expanded for the placement of 2 additional implant systems (**frontier PEAK** and **avantgard PEAK**).

Kit quirúrgico gmi frontier/avantgard - gmi frontier/avantgard surgical kit

Descripción - Description		Un.	Ref.
Fresa lanza	Lance-shaped drill	1	KYFoC1225
Fresa bola ø2,0 mm	Ball drill ø2.0 mm	1	KYFoC1502
Fresa piloto ø2,50 mm	Pilot drill ø2.50 mm	1	KYFoC2250
Fresa escalonada ø2,50 / ø3,00 mm	Step drill ø2.50 / ø3.00 mm	1	KYFoC5130
Fresa escalonada ø3,00 / ø3,50 mm	Step drill ø3.00 / ø3.50 mm	1	KYFoC5135
Fresa escalonada ø3,50 / ø4,00 mm	Step drill ø3.50 / ø4.00 mm	1	KYFoC5140
Fresa escalonada ø4,00 / ø4,50 mm	Step drill ø4.00 / ø4.50 mm	1	KYFoC5145
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4336
Fresa cortical ø4,00 mm	Cortical drill ø4.00 mm	1	KYFoC4340
Fresa cortical ø4,50 mm	Cortical drill ø4.50 mm	1	KYFoC4345
Fresa cortical ø4,90 mm	Cortical drill ø4.90 mm	1	KYFoC4349
Fresa cortical ø5,90 mm	Cortical drill ø5.90 mm	1	KYFoC4359
Llave transportador plana 45°	Manual wrench 45°	1	KYLoC0093
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave-transportador HEX-2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYLoF0096
Llave-transportador HEX-2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYLoF0111
Llave-transportador HEX-2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYLoF0091
Llave-transportador Avantgard HEX-3,00 mm IP corta	Short HEX-3.00 mm Avantgard implant carrier TI wrench	1	KYLoF0137
Llave-transportador Avantgard HEX-3,00 mm IP larga	Long HEX-3.00 mm Avantgard implant carrier TI wrench	1	KYLoF0138
Llave-transportador Avantgard HEX-3,00 mm CA	HEX-3.00 mm Avantgard implant carrier HP wrench	1	KYLoF0136
Llave carraca corta HEX-1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Llave carraca larga HEX-1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLoF0129
Alargador fresas	Drill extender	1	KYLoC0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLoC0078
Indicador profundidad ø2,50 mm	Depth gauge ø2.50 mm	1	KYLoC0165

Kit fresas gmi PEAK - gmi PEAK drill kit

Descripción - Description		Un.	Ref.
Fresa escalonada ø2,00 / ø2,80 mm	Step drill ø2.00 / ø2.80 mm	1	KYFoC5128
Fresa escalonada ø3,20 / ø3,65 mm	Step drill ø3.20 / ø3.65 mm	1	KYFoC5137
Fresa escalonada ø3,70 / ø4,15 mm	Step drill ø3.70 / ø4.15 mm	1	KYFoC5142
Fresa escalonada ø4,15 / ø4,65 mm	Step drill ø4.15 / ø4.65 mm	1	KYFoC5147
Fresa cónica ø3,20 mm	Conical drill ø3.20 mm	1	KYFoC7032
Fresa cónica ø3,70 mm	Conical drill ø3.70 mm	1	KYFoC7037
Fresa cónica ø4,20 mm	Conical drill ø4.20 mm	1	KYFoC7042
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4337
Fresa cortical ø4,15 mm	Cortical drill ø4.15 mm	1	KYFoC4342
Fresa cortical ø4,65 mm	Cortical drill ø4.65 mm	1	KYFoC4347

Llaves de inserción avantgard PEAK vendidas por separado - Avantgard PEAK insertion wrenches sold separately

Kit quirúrgico gmi PEAK - gmi PEAK surgical kit

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gmi PEAK surgical kit consist of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi frontier PEAK** and **avantgard PEAK** dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134°C.

To complement this kit, **gmi** has the **frontier/avantgard drill Kit**, allowing it to be easily expanded for the placement of 2 additional implant systems (**frontier** and **avantgard**).

KYCoF3017

KYCoF3045

Kit quirúrgico gmi PEAK - gmi PEAK surgical kit

Descripción - Description		Un.	Ref.
Fresa bola ø2,0 mm	Ball drill ø2.0 mm	1	KYFoC1502
Fresa inicial ø2,0 mm	Initial drill ø2.0 mm	1	KYFoC1401
Fresa cónica ø3,20 mm	Conical drill ø3.20 mm	1	KYFoC7032
Fresa cónica ø3,70 mm	Conical drill ø3.70 mm	1	KYFoC7037
Fresa cónica ø4,20 mm	Conical drill ø4.20 mm	1	KYFoC7042
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4337
Fresa cortical ø4,15 mm	Cortical drill ø4.15 mm	1	KYFoC4342
Fresa cortical ø4,65 mm	Cortical drill ø4.65 mm	1	KYFoC4347
Fresa escalonada ø2,00 / ø2,80 mm	Step drill ø2.00 / ø2.80 mm	1	KYFoC5128
Fresa escalonada ø3,20 / ø3,65 mm	Step drill ø3.20 / ø3.65 mm	1	KYFoC5137
Fresa escalonada ø3,70 / ø4,15 mm	Step drill ø3.70 / ø4.15 mm	1	KYFoC5142
Fresa escalonada ø4,15 / ø4,65 mm	Step drill ø4.15 / ø4.65 mm	1	KYFoC5147
Indicador profundidad ø2,00 mm	Depth gauge ø2.00 mm	1	KYLoC0164
Indicador profundidad cónico ø3,20 mm	Conical depth gauge ø3.20 mm	1	KYLoC0159
Indicador profundidad cónico ø3,70 mm	Conical depth gauge ø3.70 mm	1	KYLoC0160
Indicador profundidad cónico ø4,20 mm	Conical depth gauge ø4.20 mm	1	KYLoC0161
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave Avantgard CA corta	Short Avantgard HP wrench	1	KYLoC0151
Llave Avantgard carraca corta	Short Avantgard ratchet wrench	1	KYLoC0157
Llave Avantgard carraca larga	Long Avantgard ratchet wrench	1	KYLoC0162
Llave carraca corta HEX - 1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Llave carraca larga HEX - 1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLoF0129
Alargador fresas	Drill extender	1	KYLoC0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLoC0078
Llave-transportador HEX - 2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYLoF0096
Llave-transportador HEX - 2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYLoF0111
Llave-transportador HEX - 2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYLoF0091

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